

Product Information


BIOREBA Homogenizers and Extraction Bags

Our homogenizers are ideal for tissue extraction for samples to be used in ELISA and PCR. Samples are contained within extraction bags throughout the sample processing, thereby eliminating cross-contamination and the need of washing between samples.

HOMEX 6

Art. No. 400005 / 400004


- Semi-automated homogenizer
- Available with 230 V or 115 V external power supply
- Tissue samples of different thickness may be homogenized
- About 250 samples per hour can be processed
- Can be equipped with two different racks (see below) to hold different sizes of BIOREBA extraction bags (see verso)
- Changing the rack is fast and easy



Standard rack

Art. No. 400014


- Our HOMEX 6 is equipped with the standard rack
- Holds 1-5 BIOREBA extraction bags Universal or Universal-U-form



Rack 4/25

Art. No. 400016

- Optional rack for HOMEX 6
- Holds 1-4 BIOREBA extraction bags Universal-Long



Homogenizer hand model

Art. No. 400010

- Suitable for small sample numbers and for samples of softer tissue (leaves, sprouts, flowers etc.)
- Can be used with all types of BIOREBA extraction bags



HOMEX XS

NEW

Art. No. 400008

- Automated easy-to-use system
- Space saving and low weight machine
- High efficiency with homogenization of more than 100 samples per hour
- For homogenization of all kinds of tissues without sharp edges
- Filtered extract can be used for analyses such as ELISA, PCR or other analytical applications
- Optimized for BIOREBA extraction bags Universal
- Adjustable rotation speed
- Adjustable homogenization time
- Stop timer function
- Ergonomic handling

