


EQUINE 

DISEASE / PATHOGEN	PRODUCT NAME	METHOD	PRODUCT CODE	FORMAT	REACTIONS
Covid-19 / SARS-CoV-2	ID Screen® SARS-CoV-2 Double Antigen Multi-species	Double antigen ELISA	COVIDA-2P	2 plates	192
			COVIDA-5P	5 plates	480
Equine Infectious Anemia (EIA)	ID Screen® Equine Infectious Anemia Double Antigen	Double antigen ELISA	EIAED-1P	1 plate	96
			EIAED-4P	4 plates	384
	Positive freeze-dried horse serum		MRI-EIA	1 ml	
	IDvet EIA AGID	AGID (Coggins test)	EIA-AGID	With agar gel	200
			EIA-AGID-NOGEL	Without agar gel	200
	IDvet agar gel for EIA AGID	AGID	GEL-EIA	200 ml	up to 160
	EIA AGID CUTTER	AGID	EIACUTTER	-	
Positive freeze-dried horse serum		MRI-EIA-AGID	0.5 ml		
Equine Viral Arteritis	ID Screen® Equine Viral Arteritis Indirect	Indirect ELISA Biwell format	EVAB-4P	4 plates	192
	Positive freeze-dried horse serum		MRI-EVA	0.5 ml	
Glanders	ID Screen® Glanders Double Antigen Multi-species NEW	Double Antigen ELISA	GLANDA-5P	5 plates	480
	Freeze-dried equine inactivated serum		MRI-GLAN	0.5 ml	
Influenza A	ID Screen® Influenza A Antibody Competition Multi-species	Competitive ELISA	FLUACA-2P	2 plates	192
			FLUACA-5P	5 plates	480
			FLUACA-10P	10 plates	960
Strangles / Streptococcus equi	ID Screen® Streptococcus equi Indirect	Indirect ELISA	STREPTOS-2P	2 plates	192
	Positive freeze-dried horse serum		MRI-STREPTO	0.5 ml	
Trichinellosis	ID Screen® Trichinella Indirect Multi-species	Indirect ELISA	TRICHIS-MS-2P	2 plates	192
West Nile Virus	ID Screen® West Nile Competition Multi-species	Competitive ELISA	WNC-1P	1 plate	96
			WNC-2P	2 plates	192
	IgG-positive freeze-dried horse serum		MRI-WN	1 ml	
	ID Screen® West Nile IgM Capture	Capture ELISA Biwell format	WNIGM-1P	1 plate	48
			WNIGM-4P	4 plates	192
IgM-positive freeze-dried horse serum		MRI-WNIGM	0.5 ml		